

ASUS LU800 VET (for veterinary care)

Ultrasound Scans Anytime, Anywhere






ASUS Handheld Ultrasound for Veterinary Care

Tough, durable and reliable

Innovative ASUS handheld ultrasound technologies empower veterinarians to provide better care for patients. **These compact** devices are easy to carry and ideal for point-of-care testing in the clinic, as well as farm and house calls. From small pets to large livestock, we provides real-time high-definition imagery to allow for quick and confidence diagnosis.



Emergency procedures



Small pets



Large livestock



Product Highlights

5

Probes

2.5_{hrs}

Long Battery Life

290g

Light-weighted design

5G Wi-Fi

True wireless design

4

OS supported

5

Function modes

16

Presets

20+

Parameters

IP44/67

Water-resistant



Probes Characteristics




	Probes	LU800L	LU800C	LU800M	LU800PA	LU800E	
	Types of array	Linear	Convex	Microconvex	Phased Array	Endocavity	
	Frequency	4.2-12.5MHz	2-5MHz	3.6-8.5MHz	1.7-3.7MHz	3.6-8.5MHz	
	Max Depth	12.6 cm	30 cm	12.6 cm	30 cm	15 cm	
	Elements	128	128	128	64	128	
	Field of View	NA	60°	100°	90°	151°	
	Radius	NA	60 mm	15 mm	NA	10 mm	
	Dimensions	168x58x34	177×71×34	176x58x34	174x58x34	345x58x34	
ı	Weight	270 g	305 g	265 g	280 g	295 g	
	Image Modes	B/M/CF/PW/PD					

JPEG/PNG/BMP/MP4/DICOM

3000 mAh

Image File


Battery







Images Modes









B mode

CF mode

M mode

PW mode



Images Parameters

THI

Freq

FPS

M PRF

Color PRF

Color Wall filter

Angle

Spatial Compd.

Enhancement

Color Gain

Gate

Reverse

Gain

Gray Map

TGC

Steering Angle

PRF

Gain

Dynamic Range

Persistence

Invert

Color Threshold

Wall Filter



Al Voice Control

Hands-free access to key imaging functions

Single operator

Eliminate the need for additional assistant or second pair of hands

Efficiency Upgrade

Streamline your workflow—no more reaching for the screen.

Sterile Safety Boost

Safeguard your sterile field with voice-controlled wireless scanning, eliminating the risk of screen contact.

Independent Operation Advantage

No extra help needed; simplify tasks by eliminating the need for an additional assistant or a second pair of hands.

User scenario : Give the voice command to pad/phone/desktop/laptop → The APP will adjust the parameters automatically according to the voice command

Commands of AI voice control

Function	Voice Command		
Capture Image	Capture Image		
Capture Video	Capture Video		
Increase Depth	Increase Depth		
Decrease Depth	Decrease Depth		
Increase Gain	Increase Gain		
Decrease Gain	Decrease Gain		
Freeze	Freeze		
Center Line	Center Line		
Needle Enhance	Needle Mode		
B Mode	B Mode		
M Mode	M Mode		
CF Mode	Color Mode		
PD Mode	Power Mode		
PW Mode	Pulsed Wave Mode		


Live demo video





ASUS Ultrasound Cloud

Organize and access the exams anytime, anywhere



Ultrasound image upload

Send the exam information to ASUS Ultrasound Cloud after performing the scanning.

Access the exams anywhere

Review the exams stored in the Cloud even you're not in the hospital or clinic


Share for telecare scenario

Make notes and send exams for remote collaboration to have comprehensive diagnosis and better care to patients



Testing Standards

ltem	Information			
Battery	 UN 38.3, Lithium Battery Transportation EN IEC 62133 			
Acoustic Standards	• EN IEC 60601-2-37:2008/AMD1:2015, Medical electrical equipment, Part 2-37, Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment			
Biocompatibility Standards	 EN ISO 10993-1:2009, Biological evaluation of medical devices EN ISO 10993-5:2009, Biological evaluation of medical devices ISO 10993-10:2010, Biological evaluation of medical devices 			
Safety Standards	 IEC 60601-1:2005+AMD1:2012 / EN 60601-1:2006+ A1 2013 CSV Medical electrical equipment, Part 1: General requirements for basic safety and essential performance IEC 60601-1-2: 2014 / EN 60601-1-1:2015 Medical electrical equipment, Part 1-2: General requirements for basic safety and essential performance, Collateral Standard: Electromagnetic Capability - Requirements and tests AIUM/NEMA UD 2- 2004 2009 NEMA Standards Publication UD 2-2004 (R2009) Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment, Revision 3. (Radiology)AIUM/NEMA UD 3- 2004 2009 NEMA Standards Publication UD 3-2004 (R2009) Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment EN IEC 62304 2006 Medical device software - Software life cycle processes IEC 62366-1: 2015/EN 62366-1:2015 Medical devices, Application of usability engineering to medical devices IEC 60601-1-6 / EN 60601-1-6 Usability ISO 15223-1 2016 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied ISO 13485:2016 Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes 			





Awards of ASUS Handheld Ultrasound Solution (LU800)



"In search of good design" - the Red Dot Design Award is one of the world's largest design competitions.



The UI design of ASUS portable ultrasound solution won the iF Design Award 2023, a global symbol of excellent design.



ASUS Portable Ultrasound won the Taiwan Excellence Award, a world-renowed design award.



The product design of ASUS Portable Ultrasound Solution won the iF Design Award 2022, a global symbol of excellent design



Product Highlights

Probes
Linear/Convex

4hrs
Long Battery Life

2_M
Cable Length

5G Wi-Fi

True wireless design

4

OS supported

5

Function modes

6

Presets

20+

Parameters

IP67

Water-resistant



Handheld Ultrasound

Special version for livestock



Thank You